




TUFF-POST (Square) with Collared Eyebolts

- ▲ TUFF-POST is designed to raise the height of a height safety system using either a 600mm or 900mm TUFF-POST.
- ▲ TUFF-POST is a sealed unit made entirely from stainless steel.
- Designed to be submerged into garden beds, water and other mediums.

The building or structure for the anchorages must be assessed by an engineer, unless it is clear to a competent height safety installer that the structure is adequate. Refer to system information for site specific use.



Rated 24kN for 2x12kN anchors (two person anchor)

- 1. To comply with Australian/European Standards each TUFF-POST unit must be tested after installation.
- 2. Test consists of ultimate pull out force proof loading to 50% (12kN) of design purpose (24kN) of anchorage.
- 3. Allow sufficient time for curing before testing, see instructions on fastener container.
- 4. Recommended chemical anchor: Fischer FIS-V as per Fischer Product Supplement Data sheets.

COMPLETE UNIT - PRODUCT CODE: TUFF-POST600_24-C-2COL				
Product Parts	Description	Qty	Material	
EYEBOLT M16x26mm	Collared Eyebolt: M16x26mm	2	316SS	
TUFF-POST600_24	TUFF-POST 600mm	1	316SS	
CON-M16x200	Anchor Rod M12X200, M12 Spring Washer, M12 Nut, M12 Washer	4	316SS	



safetylink.com





Safety Centre Training

COMPLETE UNIT - PRODUCT CODE: TUFF-POST900_24-C-2COL				
Product Parts	Description	Qty	Material	
EYEBOLT M16x26mm	Collared Eyebolt: M16x26mm	2	316SS	
TUFF-POST900_24	TUFF-POST 900mm	1	316SS	
CON-M16v200	Anchor Rod M12X200 M12 Spring Washer M12 Nut M12 Washer	1	31655	


All products must be installed in accordance with SafetyLink's installation handbook, using only products supplied by SafetyLink Pty Ltd.











Patents: SafetyLink Pty Ltd has a multitude of patents, patents pending, design applications, trademarks and copyrighted documents both lodged and issued. Should you wish to know the progress of our intellectual property on a specific product please email us on ip@safetylink.com and quote the product code.

Doc ID: MKT-1374 Version: 2.0

